

# RESafety Analysis

## Insulin Infusion System

**Author:** Moniky Ribeiro

**Date:** June 2025

**Project or Institute:** CIn/UFPE

**System:** Insulin Infusion Pump (IIP)

**Iteration:** 1<sup>a</sup>

## Step 1 - Define Safety-Critical System (SCS) Scope

### 1.1. General Concerns

#### 1.1.1 Analysis Objectives

The purpose of this analysis is to model an Insulin Infusion Pump (IIP) through the iterative RESafety process, generating successive refinements of the system's safety analysis artifacts.

#### 1.1.2 System Definition

The Insulin Infusion Pump (IIP), a safety-critical system, is designed to support the treatment of Type 1 Diabetes Mellitus. Automated IIPs enhance treatment flexibility by managing multiple stages of insulin delivery, effectively mimicking physiological responses. These devices administer both rapid-acting (bolus) and continuous (basal) insulin doses.

#### 1.1.3 Resources Needed for Analysis

- Articles:
  - Martinazzo (2022);
  - Martins et al. (2015);
  - Zhang et al. (2011, 2010);
  - Bas (2020);
  - Gonzalez Atienza et al. (2024)
- Books
  - Leveson & Thomas (2018);
  - Martins & Gorschek (2021)
- General Guidelines and Manuals

#### 1.1.4 System Boundary

The system boundary encompasses activities from the moment the patient configures the infusion settings until the correct dosage is delivered via the catheter.

#### 1.1.5 Components

- Patient
- Infusion Insulin Pump
- Infusion Set

## 1.2. Safety Concerns

### 1.2.1 Identify Accidents

- **A1** - Risk of death
- **A2** - Risk of injury

### 1.2.2 Identify System-Level Hazards

- **H1** - Hypoglycemia [A1, A2]
- **H2** - Hyperglycemia [A2]

### 1.2.3 Identify System Constraints

- **SC-01** - The system must not administer insulin in excess of the prescribed dose or in unintended circumstances. [H1]
- **SC-02** - The system must ensure that the prescribed insulin dose is delivered at the correct time and in the correct amount. [H2]

### 1.2.4 Define the responsibilities

Entity	Responsability
E1 – Patient (Human Controller)	<p><b>R-01:</b> Ensure that infusion settings are correctly configured and correspond to the medical prescription [SC-01, SC-02]</p> <p><b>R-02:</b> Verify that the device interface confirms the programmed dose before administration [SC-01]</p>
E2 - Insulin Infusion Pump	<p><b>R-03:</b> Administer insulin only according to validated infusion parameters and prevent unauthorized dosages [SC-01]</p> <p><b>R-04:</b> Monitor timing and quantity of delivery to ensure correct dose is given at the right time [SC-02]</p> <p><b>R-05:</b> Detect anomalies (e.g., occlusions, over-delivery) and alert the user immediately [SC-01, SC-02]</p>


E3 - Infusion Set	<b>R-06:</b> Maintain physical integrity to prevent leaks or unintended flow of insulin [SC-01]  <b>R-07:</b> Ensure correct and timely delivery of insulin from pump to patient [SC-02]
E4 - Patient (Human Body)	<b>R-08:</b> Respond physiologically to insulin in a way that is consistent with treatment expectations (acknowledging variability) [SC-02]

### 1.2.5. Other Artifacts


- *Not applicable*

## Step 2 - Define the iStar4Safety Models


- **SD Model**



- **SR Model**



### Step 3 - Define the Control Structure



### Step 4 - Identify Unsafe Control Actions (UCAs)

Control Action	From/To	Not Providing Causes Hazard	Providing Causes Hazard	Too Early, Too Late, Out of Order	Stopped Too Soon, Applied Too Long
Program insulin dosage (R-1)	Patient / Infusion Pump	<b>UCA-01:</b> Patient does not provide "Program insulin dosage" when insulin is required, leading to underdose [H1]	<b>UCA-02:</b> Patient provides "Program insulin dosage" with a value higher than prescribed, leading to overdose [H2]  <b>UCA-03:</b> Patient provides "Program insulin dosage" with a value lower than	<b>UCA-04:</b> Patient provides "Program insulin dosage" too late, leading to hyperglycemia [H1]  <b>UCA-05:</b> Patient provides "Program insulin dosage" too early, leading to premature insulin	<i>Not applicable</i>

			prescribed, leading to underdose [H1]	administration and resulting in hypoglycemia [H2]	
--	--	--	---------------------------------------	---	--

### Additional hazards cause identified independently of the STPA results

Hazard Cause
<b>HC-01:</b> The pump is misplaced or inaccessible to the patient.[H2]

## Step 5 - Identify Controller Constraints

Unsafe Control Action	Controller Constraint
<b>UCA-01:</b> Patient does not provide “Program insulin dosage” when insulin is required, leading to underdose. [H1]	<b>C-01:</b> The patient must program the insulin dosage whenever insulin is required, according to clinical guidance. [UCA-01]
<b>UCA-02:</b> Patient provides “Program insulin dosage” with a value higher than prescribed, leading to overdose. [H2]	<b>C-02:</b> The patient must ensure the programmed insulin dosage does not exceed the value prescribed by the physician. [UCA-02]
<b>UCA-03:</b> Patient provides “Program insulin dosage” with a value lower than prescribed, leading to underdose. [H1]	<b>C-03:</b> The patient must verify that the programmed dosage meets the minimum prescribed threshold to avoid underdosing. [UCA-03]
<b>UCA-04:</b> Patient provides “Program insulin dosage” too late, leading to hyperglycemia. [H1]	<b>C-04:</b> The patient must program the insulin dosage in a timely manner, according to the prescribed administration window. [UCA-04]
<b>UCA-05:</b> Patient provides “Program insulin dosage” too early, leading to premature insulin administration and resulting in hypoglycemia. [H2]	<b>C-05:</b> The patient must not program the insulin dosage before the appropriate physiological or dietary conditions occur. [UCA-05]
<b>HC-01:</b> The pump is misplaced or inaccessible to the patient.	<b>C-06:</b> The insulin pump must always be correctly placed and readily accessible to the patient.

## Step 6 - Analyze Loss Scenarios and Derive Safety Requirements

UCA	Loss Scenario (LS)	Safety Requirement (SR)
-----	--------------------	-------------------------

<p><b>UCA-01:</b> Patient does not provide “Program insulin dosage” when insulin is required, leading to underdose [H1]</p>	<p><b>LS-01:</b> The patient forgets to program the dose after the meal, resulting in hyperglycemia. [UCA-01] <i>Martinazzo (2022)</i></p> <p><b>LS-02:</b> The system does not issue a reminder to program the dose after detecting a meal event. [UCA-01] <i>Ribeiro et al. (2024)</i></p>	<p><b>SR-01:</b> The system shall generate an alert if insulin is not programmed within 15 minutes after a meal is detected. [LS-01] <i>Zhang et al. (2011)</i></p> <p><b>SR-02:</b> The interface must maintain a visible warning if no insulin programming is detected post-meal. [LS-02] <i>Ribeiro et al. (2024)</i></p>
<p><b>UCA-02:</b> Patient provides “Program insulin dosage” with a value higher than prescribed, leading to overdose [H2]</p>	<p><b>LS-03:</b> The patient repeats a bolus due to lack of feedback on recent insulin administration. [UCA-02] <i>Zhang et al. (2010)</i></p> <p><b>LS-04:</b> Patient misinterprets the prescribed dose and enters a value higher than medically indicated. [UCA-02] <i>Zhang et al. (2011)</i></p>	<p><b>SR-03:</b> The system shall display recent insulin activity clearly before accepting a new dose. [LS-03] <i>Zhang et al. (2011)</i></p> <p><b>SR-04:</b> The system shall cross-check manual input with prescription data and alert if excess dosage is detected. [LS-04] <i>Zhang et al. (2011)</i></p>
<p><b>UCA-03:</b> Patient provides “Program insulin dosage” with a value lower than prescribed, leading to underdose [H1]</p>	<p><b>LS-05:</b> The patient reduces the dose to avoid hypoglycemia without clinical basis. [UCA-03] <i>Martinazzo (2022)</i></p> <p><b>LS-06:</b> The system does not notify that the entered dose is below clinical expectation. [UCA-03] <i>Zhang et al. (2011)</i></p>	<p><b>SR-05:</b> The system must recommend confirmation when the user’s dose is significantly below the recommended amount. [LS-05, LS-06] <i>Zhang et al. (2011)</i></p>
<p><b>UCA-04:</b> Patient provides “Program insulin dosage” too late, leading to hyperglycemia [H1]</p>	<p><b>LS-07:</b> The patient delays programming due to being busy or distracted, compromising glycemic control. [UCA-04] <i>Martinazzo (2022)</i></p> <p><b>LS-08:</b> The system accepts bolus entry even after blood glucose spike already occurred. [UCA-04] <i>Ribeiro et al. (2024)</i></p>	<p><b>SR-06:</b> The interface must issue periodic prompts for pending bolus if blood glucose remains elevated and no dose is scheduled. [LS-07] <i>Zhang et al. (2011)</i></p> <p><b>SR-07:</b> The system must block bolus entries considered ineffective post-prandial, requiring physician override. [LS-08] <i>Ribeiro et al. (2024)</i></p>

<p><b>UCA-05:</b> Patient provides "Program insulin dosage" too early, leading to premature insulin administration and resulting in hypoglycemia [H2]</p>	<p><b>LS-09:</b> Patient programs insulin and forgets to eat, leading to insulin drop without carbohydrate intake. [UCA-05] <i>Zhang et al. (2010)</i></p> <p><b>LS-10:</b> Patient assumes a meal is imminent, but it is delayed due to unforeseen events. [UCA-05] <i>Martins et al. (2015)</i></p>	<p><b>SR-08:</b> System must require user confirmation that the meal is occurring before completing bolus delivery. [LS-09] <i>Zhang et al. (2011)</i></p> <p><b>SR-09:</b> If a meal confirmation is not received within a set time, bolus must be suspended or canceled automatically. [LS-10] <i>Martins et al. (2015)</i></p>
---	---	---

### Additional hazards cause identified independently of the STPA results

Hazard Cause	Loss Scenario	Safety requirement
<p><b>HC-01:</b> The pump is misplaced or inaccessible to the patient.</p>	<p><b>LS-11:</b> The patient is in a critical condition and does not remember where the pump was placed.</p>	<p><b>SR-10:</b> The pump must have an associated mobile application that allows a "locate pump" function to trigger an audible alarm when activated.</p>



## Step 7 - Update the iStar4Safety Models

